

WESTERN TRADING COMPANY

Absorbable/Non-Absorbable Fixation Devices



Proven to perform

Our fixation portfolio is comprehensive in scope and singular in reputation. These market leading products are designed to support durable hernia repairs while offering the confidence that comes only with vast clinical experience. From absorbable fixation, to permanent tacks and staples in a wide variety of sizes and loading options, this versatile family of fixation devices answers the specific needs of surgeons and is designed with the intent to optimize clinical performance.

RELIATACK™ FIXATION DEVICE

The ReliaTack™ fixation device offers one-of-a-kind articulating technology, interchangeable reloads, and screw-like tack designs for hernia repair.

ABSORBATACK™ FIXATION DEVICE

The AbsorbaTack™ 5 mm fixation device is a sterile, single use device for fixation of prosthetic material, such as hernia mesh, to soft tissue. The tack is constructed of an absorbable synthetic polyester copolymer derived from lactic and glycolic acid.

PROTACK™ FIXATION DEVICE

The ProTack™ 5 mm fixation device is a sterile, single use device for fixation of prosthetic material, such as hernia mesh, to soft tissue. The tack is helical and made of titanium. The ProTack™ device comes in a 30 tack configuration.

ENDO UNIVERSAL™ STAPLER

The Endo Universal™ stapler has application in endoscopic hernia repair surgery procedures for fixation of mesh. The device can be rotated 360° around its axis and can be used in four different positions: unarticulated position at 0°, and articulated 20°, 45°, and 65°.

MULTIFIRE ENDO HERNIA™ STAPLER

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.

The MultiFire Endo Hernia™ stapler is designed for introduction and use through all appropriately sized trocar sleeves or larger sized trocar sleeves with the use of a converter. The size of the staples is determined by the selection of the 4mm or the 4.8mm single use loading kit (SULU).

STAT TACK™ FIXATION DEVICE

The Stat Tack™ 5 mm fixation device delivers helical titanium tacks for fixation of mesh in open hernia repair.

TACKER™ FIXATION DEVICE

The Tacker™ 5 mm fixation device delivers helical titanium tacks for laparoscopic fixation of mesh in hernia repair.

MULTIFIRE VERSATAACK™ STAPLER

The MultiFire VersaTack™ stapler has application in general surgery procedures for fixation of mesh in the repair of hernia defects and in other surgical specialties for approximation of tissue(s).

At a glance

Brand / Division	Western Trading Company
Specialties	GI Surgery
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad



Scan to view online



Scan for WhatsApp enquiry