

WESTERN TRADING COMPANY

Absorbable Sutures



Our absorbable sutures provide strength retention and absorption profiles to meet your surgical specialty and tissue healing requirements.

CAPROSYN™ MONOFILAMENT ABSORBABLE SUTURES

Caprosyn™ synthetic absorbable suture's quick absorption (56 days), durable polymer material and monofilament construction provide smooth passage through tissue, better approximation or ligation, and minimized potentiation of infection. Patients can experience less discomfort during healing, less wound complications and a faster recovery.

POLYSORB™ BRAIDED ABSORBABLE SUTURES

Polysorb™ sutures are composed of Lactomer™ glycolide/lactide copolymer, which is a synthetic polyester composed of glycolide and lactide (derived from glycolic and lactic acids). They are prepared by coating the suture with a mixture of a caprolactone/glycolide copolymer and calcium stearoyl lactylate. The sutures are colored violet to increase visibility and are also available undyed.

BIOSYN™ MONOFILAMENT ABSORBABLE SUTURES

Biosyn™ synthetic absorbable sutures are prepared from Gylcomer™ 631, a synthetic polyester composed of glycolide (60%), dioxanone (14%), and trimethylene carbonate (26%). These sutures are colored violet to increase visibility and are also available undyed.

MAXON™ MONOFILAMENT ABSORBABLE SUTURES

Maxon™ synthetic absorbable sutures are prepared from polyglyconate, a copolymer of glycolic acid and trimethylene carbonate.

MAXON™ CV MONOFILAMENT ABSORBABLE SUTURES

Maxon™ CV synthetic absorbable sutures are prepared from polyglyconate, a copolymer of glycolic acid and trimethylene carbonate.

At a glance

Brand / Division	Western Trading Company
Specialties	Oncology, Bariatric Surgery, Urology & Transplant, Neuro & Spine
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESCULAP EinsteinVision 3.0 FI · Needle Holders · Bi-Polar Lap Instruments · LigaSure™ Maryland Jaw Sealer-Divider

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.