

MEDTRONIC ENDOSCOPY

Barrx™ Radiofrequency Ablation System



The Barrx™ radiofrequency ablation system treats patients with Barrett's esophagus by eradicating diseased epithelium with precisely controlled radiofrequency energy. Built around the Barrx™ flex RFA energy generator, the platform offers a full range of balloon and focal catheters for both circumferential segments and small focal areas.

- Barrx™ flex RFA energy generator delivers controlled bipolar radiofrequency energy
- Barrx™ 360 Express RFA balloon catheter for circumferential ablation of larger treatment areas
- Barrx™ 60 and 90 RFA focal catheters for short and intermediate segments or touch-up
- Barrx™ channel RFA endoscopic catheter delivers ablation through the working channel of a flexible endoscope
- Clinically proven therapy for eradication of Barrett's epithelium

At a glance

Brand / Division	Medtronic Endoscopy
Specialties	Endoscopy
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

GI Genius™ Intelligent Endoscopy Module · PillCam™ SB 3 Capsule Endoscopy System · PillCam™ Crohn's System · Bravo™ Calibration-Free Reflux Testing System · Beacon™ EUS Delivery System · Endoflip™ 300 Impedance Planimetry System

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.