

MEDTRONIC ENDOSCOPY

## Beacon™ EUS Delivery System



The Beacon™ endoscopic ultrasound delivery system is a tissue acquisition platform featuring interchangeable needles — including Beacon™ FNA, Beacon™ FNF fiducial and SharkCore™ FNB needles — designed for consistent, predictable acquisition of tissue samples with intact cellular architecture.

- Interchangeable FNA, FNB and fiducial needles on one delivery system
- Exchange needles without removing the delivery system, improving procedural workflow
- SharkCore™ FNB needle design supports histology-quality cores with intact cellular architecture
- Built-in safety features help protect staff from needle-stick injury
- Supports EUS-guided sampling of pancreatic and other GI lesions

### At a glance

Brand / Division	Medtronic Endoscopy
Specialties	Endoscopy
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

### Also available in this specialty

GI Genius™ Intelligent Endoscopy Module · PillCam™ SB 3 Capsule Endoscopy System · PillCam™ Crohn's System · Bravo™ Calibration-Free Reflux Testing System · Barrx™ Radiofrequency Ablation System · Endoflip™ 300 Impedance Planimetry System

### Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

## Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.