

AESCULAP

Maryland AdTec Monopolar



For years, Aesculap® has been actively involved in the development of laparoscopic techniques in the instruments sector, introducing products such as AdTec® single use, AdTec® needle holders as well as the Multi-Fire Clip Applier Challenger® Ti-P.

Features

- “Universal”: modular design consisting of four components, facilitates the disassembling, reprocessing, maintenance and sterile supply.
- “One-Click” design: assembly and disassembly in one step only.
- “PEEK Isolation”: Enhanced spark protection and resisting characteristic. Wear-resistant and sterilization-proof.
- “Ratchet mechanism”: ratchet can be deactivated.
- “Jaw mechanism”: non-protruding joints during opening and closure of jaw parts and thus a lower risk that tissue or suture material gets caught in the jaw mechanism.
- “Cleaning/ sterilisation”: complete disassembly of the instrument into 4 Parts leading to better cleaning of the instrument.

At a glance

Brand / Division	Aesculap
Specialties	Gynecology, Urology & Transplant
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESCULAP EinsteinVision 3.0 FI · “Non-slip” Lap Needle Holders · Bi-Polar Lap Instruments · LigaSure™ Maryland Jaw Sealer-Divider

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.