


MEDTRONIC ENDOSCOPY

## Bravo™ Calibration-Free Reflux Testing System



The Bravo™ calibration-free reflux testing system is a catheter-free, capsule-based pH monitoring system used to determine objectively whether acid reflux is causing a patient's symptoms. A small capsule attached to the esophageal lining wirelessly transmits pH data to a worn recorder, allowing patients to maintain normal diet and activities during the test.

- Up to 96 hours of reflux monitoring — the longest available wireless pH test window
- Calibration-free capsule simplifies setup and workflow
- Capsule can be placed after either endoscopy or manometry
- Patient-worn recorder captures reflux events, symptoms and activity data
- Catheter-free design improves patient comfort versus traditional transnasal pH testing

### At a glance

Brand / Division	Medtronic Endoscopy
Specialties	Endoscopy
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

### Also available in this specialty

GI Genius™ Intelligent Endoscopy Module · PillCam™ SB 3 Capsule Endoscopy System · PillCam™ Crohn's System · Barrx™ Radiofrequency Ablation System · Beacon™ EUS Delivery System · Endoflip™ 300 Impedance Planimetry System

### Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

## Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.