

## VALLEYLAB™ · ELECTROSURGICAL ENERGY

# Electrosurgical Accessories

A complete range of Valleylab™ electrosurgical accessories for the FX8 and FT10 energy platforms — from value disposable pencils to the top-of-the-line Force TriVerse™, plus patient return electrodes and a full set of interchangeable electrode tips.

## Electrosurgical Pencils

	VL2610	E2516	FT3000
Tier	Value	Standard	Top of the line
Switch	Hand control, replaceable blade	Push-button (Cut / Coag)	Three-button + sterile-field slider
Energy modes	Cut & Coag	Cut & Coag	Cut, Coag & Valleylab™ mode
Blade electrode	Standard stainless steel	Stainless steel, hex-locking	Non-stick Edge™-coated, hex-locking
Sterile-field power control	—	—	Yes — power slider on the handle
Designed for	Standard Valleylab™ generators	Standard Valleylab™ generators	Optimised for FT10 & FX8
Best for	Routine, cost-sensitive procedures	Everyday monopolar cut & coag	Advanced cases needing less char & in-field control

- FT3000 (Force TriVerse™) is the top-of-the-line pencil — three modes, sterile-field power slider and a non-stick Edge™-coated blade, optimised for the FX8 & FT10 platforms.

## Patient Return Electrode

**E7507 — REM PolyHesive™ II Patient Return Electrode.** The adult-size patient return electrode (grounding pad) that safely returns electrosurgical current to the generator. It supports Valleylab's REM™ contact-quality monitoring, which automatically deactivates the generator if pad contact becomes unsafe — helping prevent patient burns.

## Electrode Tips

- **Blade electrode** (E1450G (Edge™) / E1551X (steel)) — General-purpose cutting and coagulation — the everyday workhorse tip.
- **Needle electrode** (E1452 (Edge™ coated)) — Fine, precise incisions and detailed dissection.
- **Insulated blade** (E1455 (PTFE insulated)) — Concentrates energy at the tip and protects adjacent tissue.
- **Ball electrode** (E1000 / E1002 (reusable)) — Surface coagulation and fulguration over a broader area.
- **Loop electrode** (E1005 (tungsten, reusable)) — Shaving or excising thin slices of tissue.
- **Electrode extensions** (E1502 / E1504) — Added reach for deep-cavity access.

## Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.