

AESFULAP

AESFULAP ELAN 4 Power System



The AESFULAP ELAN 4 is the current high-speed electric power system for neuro, spine, ENT and maxillofacial surgery — the successor to the older Hi-Line XS — intended for high-speed cutting, sawing and drilling of bone. Its defining feature is a powerful direct drive motor in the handpiece that transfers power directly to the burr or blade, reducing friction, overheating and chatter.

- Direct drive motor in the handpiece for direct power transfer and minimal friction
- Lightweight, ergonomic, well-balanced handpieces for comfort during long procedures
- Universal burr lengths with clear ring coding for foolproof handpiece matching
- Fully submersible handpieces for manual and mechanical cleaning
- Plug-and-play design for quick integration into existing OR workflows

At a glance

Brand / Division	Aesculap
Specialties	Neuro & Spine
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FX8 Energy Platform · Neuro Clips Instrument · AESFULAP SQ.line KERRISON Bone Punches · Leyla Retractors · Neuro Endoscopy · Bipolar Forceps Instruments

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.