

MEDTRONIC ENT

Hydrodebrider™ Endoscopic Sinus Irrigation System



The Hydrodebrider™ system is a powered endoscopic sinus irrigation system for treating chronic sinusitis, used in the operating room or office. Winner of a Medical Device Excellence Award, it delivers a rotating spray of pressurised saline shown in laboratory studies to dislodge bacteria far more effectively than traditional nasal lavage.

- Standard and frontal handpieces deliver a rotating spray of pressurised saline at 5 mL/sec
- Flex-Tip on the standard handpiece combines suction with 220° articulation
- Standard handpiece accesses the maxillary, ethmoid and sphenoid sinuses
- Slim 2.2 mm frontal handpiece reaches hard-to-access frontal sinuses
- Powered irrigation removes bacteria associated with recurrent infection in chronic rhinosinusitis

At a glance

Brand / Division	Medtronic ENT
Specialties	ENT
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this speciality

NIM Vital™ Nerve Monitoring System · StealthStation™ ENT Navigation System · StealthStation FlexENT™ Navigation System · Straightshot™ M5 Microdebrider · NuVent™ EM Balloon Sinus Dilation System · PROPEL™ Sinus Implant

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.