

WESTERN TRADING COMPANY

Leyla Retractors



The position of the flexible holding arm and the LEYLA retractor can be adjusted in three dimensions. A locking mechanism allows fixating the devices in any position.

Product Features – Flexible Holding Arm

- Stable holding arm with single-hand locking mechanism
- For endoscopes and MINOP® system, ventriculocopy or endoscope-assisted neurosurgery

Product Features – Leyla Retractors

- Maximum ease of handling
- Adjustability of clamping and holding force
- For spatula fixation
- Various clamping elements available

At a glance

Brand / Division	Western Trading Company
Specialties	Neuro & Spine
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FX8 Energy Platform · AESCULAP ELAN 4 Power System · Neuro Clips Instrument · AESCULAP SQ.line KERRISON Bone Punches · Neuro Endoscopy · Bipolar Forceps Instruments

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.