

AESFULAP

## Needle Holders



With these needle holders AESFULAP® offers a range of instruments specially developed for suturing and knotting. The axial handle allows the instruments to be rotated very easily. The integrated overload protection reduces damage to the needle holder and needle.

- Specific grip on curved and straight needles; grip on the suture thread.
- Integrated overload protection allows holding course needles as reliably as fine ones.
- Ergonomic working posture through axial handle design.
- The axial arrangement allows holding the instrument in any position.
- Simple lock with just two fixed positions following the push ballpoint principle. The needle can thus be gripped with the lock either engaged or released.
- Facilitated suturing and knotting through mutually corresponding instrument combinations.
- Versatile needle holders available in different diameters and lengths.
- All DUROGRIP needle holders with tungsten carbide inserts in jaw parts.

### At a glance


Brand / Division	Aesculap
Specialties	Oncology, Bariatric Surgery
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

### Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESFULAP EinsteinVision 3.0 FI · Bi-Polar Lap Instruments · LigaSure™ Maryland Jaw Sealer-Divider · Sonicision™ 7 Curved Jaw Cordless Ultrasonic Dissector

### Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



## Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.