

AESCULAP

“Non-slip” Lap Needle Holders



Looking for a non-slip gripping Needle Holder for a smooth suturing procedure?

With these needle holders AESCULAP® offers a range of instruments specially developed for suturing and knotting. The axial handle allows the instruments to be rotated very easily. The integrated overload protection reduces damage to the needle holder and needle.

- Specific grip on curved & straight needles; excellent grip on the suture thread.
- Integrated overload protection allows holding force of the jaw parts constant in Lap surgery.
- Ergonomic working posture through axial handle design.
- The axial arrangement allows holding the instrument in any position.
- Simple lock with just two fixed positions following the push ballpoint principle. The needle can thus be gripped with the lock either engaged or released.
- Versatile needle holders available in different lengths 31, 37 & 42 cms.
- Facilitated suturing and knotting through mutually corresponding instrument combinations.
- All DUROGRIP needle holders with tungsten carbide inserts in jaw parts.

At a glance

Brand / Division	Aesculap
Specialties	Gynecology
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESCULAP EinsteinVision 3.0 FI · Maryland AdTec Monopolar · Bi-Polar Lap Instruments · LigaSure™ Maryland Jaw Sealer-Divider

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.