

MEDTRONIC ENDOSCOPY

Nexpowder™ Endoscopic Hemostasis System



The Nexpowder™ endoscopic hemostasis system delivers a noncontact, nonthermal hemostatic adhesive powder through a catheter for the management of active and refractory nonvariceal upper GI bleeding. On contact with the bleeding site, the powder immediately forms a durable muco-adhesive gel that naturally degrades within one to three days.

- Noncontact, nonthermal and nontraumatic hemostasis for nonvariceal upper GI bleeding
- Proprietary powder-coating technology minimises catheter clogging and maintains visibility
- Self-contained delivery device requires no CO2 or air compressors
- Distinctive blue colour confirms full coverage of the lesion
- Demonstrated 94% immediate hemostasis rate with a 3.7% rebleeding rate in clinical evaluation

At a glance

Brand / Division	Medtronic Endoscopy
Specialties	Endoscopy
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

GI Genius™ Intelligent Endoscopy Module · PillCam™ SB 3 Capsule Endoscopy System · PillCam™ Crohn's System · Bravo™ Calibration-Free Reflux Testing System · Barrx™ Radiofrequency Ablation System · Beacon™ EUS Delivery System

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.