

MEDTRONIC ENDOSCOPY

PillCam™ Crohn's System



The PillCam™ Crohn's system enables direct visualisation of both the small bowel and colonic mucosa in a single, minimally invasive procedure — the only capsule platform to do so. It is designed for assessment and ongoing monitoring of mucosal lesions that may indicate Crohn's disease.

- Dual-camera capsule visualises the small bowel and colon in one noninvasive procedure
- Adaptive frame rate increases image capture from 4 up to 35 frames per second
- Indicated for adults and children from eight years of age
- Supported by the PillCam™ recorder 3, sensor belt and sensor array ecosystem
- PillCam™ patency capsule available to verify GI tract patency in patients with suspected strictures

At a glance

Brand / Division	Medtronic Endoscopy
Specialties	Endoscopy
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

GI Genius™ Intelligent Endoscopy Module · PillCam™ SB 3 Capsule Endoscopy System · Bravo™ Calibration-Free Reflux Testing System · Barrx™ Radiofrequency Ablation System · Beacon™ EUS Delivery System · Endoflip™ 300 Impedance Planimetry System

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.