


MEDTRONIC ENT

## PROPEL™ Sinus Implant



PROPEL™ sinus implants are bioabsorbable, steroid-releasing implants placed following sinus surgery to maintain patency and locally deliver mometasone furoate to the sinus mucosa. The self-expanding implant props the sinus open while eluting corticosteroid directly where it is needed, then dissolves on its own as the sinus heals.

- Two-in-one mechanism: self-expanding support plus localised drug delivery
- Delivers 370 mcg of mometasone furoate over 30 days, with drug present in mucosal tissue up to 60 days
- Structural backbone bio-absorbs over approximately 30–45 days after placement
- Non-obstructive design allows mucociliary clearance and delivery of topical rinses
- Three variants: PROPEL™ (ethmoid), PROPEL™ Mini (ethmoid/frontal opening) and PROPEL™ Contour (frontal ostia)

### At a glance

Brand / Division	Medtronic ENT
Specialties	ENT
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

### Also available in this specialty

NIM Vital™ Nerve Monitoring System · StealthStation™ ENT Navigation System · StealthStation FlexENT™ Navigation System · Straightshot™ M5 Microdebrider · NuVent™ EM Balloon Sinus Dilation System · SINUVA™ Sinus Implant

### Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

## Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.