

MEDTRONIC ENT

Straightshot™ M5 Microdebrider



The Straightshot™ M5 microdebrider is a powered surgical handpiece used to incise or remove soft tissue, hard tissue and bone during ENT procedures. Driven by Medtronic's Integrated Power Console (IPC™), it aids operative efficiency by simultaneously cutting and suctioning fluid and debris from the surgical site.

- Up to 30,000 rpm cutting speed for efficient removal of soft tissue, hard tissue and bone
- Integrated suction keeps the surgical field clear without instrument exchanges
- Precision finger wheel for one-handed rotation of the blade tip
- Sculpted, ergonomic design that accommodates multiple hand positions and sizes
- Compatible with more than 125 application-specific blades and burs, including navigated and ApexCut™ blades

At a glance

Brand / Division	Medtronic ENT
Specialties	ENT
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

NIM Vital™ Nerve Monitoring System · StealthStation™ ENT Navigation System · StealthStation FlexENT™ Navigation System · NuVent™ EM Balloon Sinus Dilation System · PROPEL™ Sinus Implant · SINUVA™ Sinus Implant

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.