

WESTERN TRADING COMPANY

Non-Absorbable Sutures



Our Nonabsorbable sutures can be used in a range of procedures including cardiovascular, ophthalmic, and neurological surgery.

DERMALON™ MONOFILAMENT NYLON SUTURES

Dermalon™ monofilament nylon sutures are inert, non-absorbable, sterile surgical sutures, which are available dyed blue to provide visibility compared to other nylon sutures.

MONOSOF™ MONOFILAMENT NYLON SUTURES

Monosof™ monofilament nylon sutures are inert, non-absorbable, sterile surgical sutures, which are available either dyed black, with Logwood extract, or undyed (clear).

SURGIPRO™ MONOFILAMENT SUTURES

Uncoated Surgipro™ sutures (clear or pigmented) are inert, non-absorbable, sterile sutures composed of an isotactic crystalline stereoisomer of polypropylene (a synthetic linear polyolefin) and polyethylene.

VASCUFIL™ MONOFILAMENT SUTURES

Vascufil™ coated, non-absorbable, monofilament surgical sutures are composed of polybutester fiber, which gives the suture

SOFSILK™ SILK SUTURES

Sofsil™ sutures are non-absorbable, sterile, non-mutagenic surgical sutures composed of natural proteinaceous silk fibers called fibroin. The braided sutures are available coated uniformly with either silicone or a special wax mixture to reduce capillarity and to increase surface lubricity, which enhances handling characteristics, ease of passage through tissue, and knot run-down properties.

SURGIPRO™ II MONOFILAMENT SUTURES

Uncoated Surgipro™ II sutures (clear or pigmented) are inert, non-absorbable, sterile sutures composed of an isotactic, crystalline stereoisomer of polypropylene and contain polyethylene.

At a glance

Brand / Division	Western Trading Company
Specialties	GI Surgery, Oncology, Bariatric Surgery, Urology & Transplant, Neuro & Spine
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESCULAP EinsteinVision 3.0 FI · Laparoscopy/Open Staplers · Signia™ Multi-Use Loading Units + Tri-Staple™ 2.0 Cartridges · Sutures

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.