


WESTERN TRADING COMPANY

V-LOC Barbed Sutures



The V-Loc™ wound closure device is a revolutionary technology that eliminates the need to tie knots, so you can close incisions up to 50% faster without compromising strength and security.^{1,2}

Its dual-angle barbs allow suture to pass smoothly through tissue in one direction — but prevent it from slipping back through soft tissue.

The barbs hold wound surfaces together and distribute tension along the entire length of the suture, without needing knots to secure the wound.

Designed for secure, fast, and effective^{3,4} incision closure for your patients.

At a glance

Brand / Division	Western Trading Company
Specialties	GI Surgery, Oncology, Bariatric Surgery, Urology & Transplant
Local support	On-site demonstration · Clinical training · Installation · After-sales service
Availability	Hospitals across Gujarat, via Western Trading Company, Ahmedabad

Also available in this specialty

Valleylab™ FT10 Energy Platform · Valleylab™ FX8 Energy Platform · AESCULAP EinsteinVision 3.0 FI · Laparoscopy/Open Staplers · Absorbable/Non-Absorbable Fixation Devices · Laparoscopy Instruments

Why Western Trading Company

- 30+ years supplying Gujarat's operating theatres — established 1995
- Trusted by 500+ hospitals including Apollo, Zydus, Sterling, HCG and CIMS
- Authorised channel partnerships with Medtronic, Aesculap B. Braun and Spatz — genuine products, local stock and responsive service



Scan to view online



Scan for WhatsApp enquiry

Western Trading Company

611-Iscon Elegance, Near Prahladnagar Garden, S.G. Highway, Ahmedabad - 380015, Gujarat, India

+91 95372 44565 · +91 98250 06579 · +91 63517 81672

contact@westerntrading.in · www.westerntrading.in

Authorised distributor spec sheet. Product names & trademarks belong to their respective manufacturers. For full manufacturer documentation, please contact us.